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
Dimensions in [cm]

INDEX OF SYMBOLS

	CE	Delivery point (only the cable, over length 5m) of electrical power supply line. 1PH+N+PE 240 Vac - 50/60Hz - 2.2kW
	UL	Delivery point (only the cable, over length 5m) of electrical power supply line. 1PH+N+GND 120 Vac - 60Hz - 2.2kW
	CE	N.4 wall mount SCHUKO socket 1PH+N+PE - 16A 240 Vac 50/60Hz - Full load current = 16A
	UL	N.4 wall mount NEMA socket 1PH+N+GND - 20A 120 Vac 60Hz - Full load current = 20A
	UL-CE	Ethernet socket RJ45 10/100/1000 Mbit connected to LAN (Local Area Network) with cable CAT. 5e
	UL-CE	Compressed Air - ISO 8573.1 class 1:2:1 recommended (but not mandatory). Pressure 6 bar; Flow rate 40 NL/min. The line must be equipped with a reachable ball valve. Female G1/4"
	UL-CE	Ø5x110cm Exhaust duct -50Pa refer to Lab - Exhaust flow rate 50m ³ /h From lab to the roof of the building

NOTES:

- The power supply cable has to be sized according to the customer's electrical distribution system feature and the rules in force to the country
- Differential protection relay needed in power supply
- Heat dissipation: 350 kcal/h

DESCRIPTION

**Manipulation Hot Cell
Mod. FLEX 2T**


TEMA SINERGIE S.p.A.
Via Malpighi, 120
48018 Faenza (RA) Italy
T. +39 0548 822893
F. +39 0548 821640
www.temasinergie.com

DRAWING CODE	SLST116 - 01	REV	00
SCALE	1:40 - 1:25		
DATE	04/10/2018	APPROVED BY	
DRAFTSMAN	M.M.	SIGN	
AS BUILT			

STANDARD CONFIGURATION
All pictures and specifications included in this document are indicative and may be subject to change without notice


CUSTOMER

FORMAT A3

SECT. Y-Y

FLOOR ARRIVAL OF RADIOACTIVE LINE




FRONT CELL - PHARMA WALL

CHANGES

DATE

DRAFTSMAN

APPROVED


	Box 1	Box door opening 2
W	92 cm	72 cm
D	82 cm	/
H	66 cm	43 cm


Dimensions in [cm]

DATE	CHANGES	DRAFTSMAN	APPROVED	CUSTOMER	FORMAT
				<p>TEMA SINERGIE</p> <p>Tema Sinergie S.p.A. Via Malpighi, 120 48018 Faenza (RA) Italy T. +39 0546 822693 F. +39 0546 821640 www.temasinergie.com</p>	A3
				<p>STANDARD CONFIGURATION</p> <p>All pictures and specifications included in this document are indicative and may be subject to change without notice</p>	
				<p>DESCRIPTION Manipulation Hot Cell Mod. FLEX 2T</p>	
				<p>DRAWING CODE: SLST116 - 02</p>	
				<p>SCALE: 1:20</p>	
				<p>DATE: 04/10/2018</p>	
				<p>DRAFTSMAN: M.M.</p>	
				<p>APPROVED BY: _____</p>	
				<p>AS BUILT</p>	
				<p>DATE: _____</p>	
				<p>SIGN: _____</p>	
				<p>Unauthorized copy usage or distribution of this drawing is forbidden under the law</p>	




44

Vial Extraction System
RIGHT SIDE

Vial Extraction System
FRONT SIDE





Dimensions in [cm]

DESCRIPTION Manipulation Hot Cell Mod. FLEX 2T				Tema Sinergie S.p.A. Via Malpighi, 120 48018 Faenza (RA) Italy T. +39 0546 822093 F. +39 0546 821640 www.temasinergie.com	
DRAWING CODE	SLST116 - 03	REV	00	CUSTOMER	
SCALE	1:40			STANDARD CONFIGURATION All pictures and specifications included in this document are indicative and may be subject to change without notice	
DATE	04/10/2018	APPROVED BY		FORMAT A3	
DRAFTSMAN	M.M.	SIGN		All pictures and specifications included in this document are indicative and may be subject to change without notice	
AS BUILT		DATE		Unauthorised copy, usage or distribution of this drawing is forbidden under the law	
DATE		DRAFTSMAN	APPROVED		
CHANGES					

62

See page 1
FLEX Installation requirement

PRECHAMBER	DESCRIPTION	TOTAL WEIGHT	FRAME AREA
		400 kg	0,43 kg

Dimensions in [cm]

INDEX OF SYMBOLS

(S1)	UL-CE	ØEXT10cm Exhaust duct -50Pa refer to Lab - Exhaust flow rate 30m³/h From lab to the roof of the building
-------------	-------	--

NOTES:
 - The power supply cable has to be sized according to the customer's electrical distribution system feature and the rules in force to the country
 - Differential protection relay needed in power supply
 - Heat dissipation: 580 kcal/h

DESCRIPTION

**Manipulation Hot Cell
Mod. FLEX 2T + PRECHAMBER**

TEMA SINERGIE S.p.A.
Via Malpighi, 120
48018 Faenza (RA) Italy
T. +39 0540 822063
F. +39 0540 821640
www.temasinergie.com

DRAWING CODE	SLST116 - 04	REV	00
SCALE	1:40 - 1:25		
DATE	04/10/2018		
DRAFTSMAN	M.M.	APPROVED BY	
AS BUILT			
DATE		SIGN	

STANDARD CONFIGURATION
 All pictures and specifications included in this document are indicative and may be subject to change without notice

CUSTOMER

Unauthorized copy usage or distribution of this drawing is forbidden under the law

DATE

CHANGES

DRAFTSMAN




APPROVED

SIGN

DATE

FORMAT

A3

	Box 1	Door opening 2	Hatch 3
W	49.4 cm	41 cm	17 cm
D	76 cm	/	/
H	69.6 cm	59.1 cm	30.3 cm

Dimensions in [cm]


DESCRIPTION Manipulation Hot Cell Mod. FLEX 2T + PRECHAMBER			
DRAWING CODE SLST116 - 05		CUSTOMER	
SCALE 1:10	REV 00	STANDARD CONFIGURATION All pictures and specifications included in this document are indicative and may be subject to change without notice	
DATE 04/10/2018			
DRAFTSMAN M.M.	APPROVED BY		
DATE AS BUILT	SIGN	All pictures and specifications included in this document are indicative and may be subject to change without notice	
DRAFTSMAN	APPROVED	FORMAT A3	
CHANGES	Unauthorized copy, usage or distribution of this drawing is forbidden under the law		

Tema Sinergie S.p.A.
 Via Malpighi, 120
 48018 Faenza (RA) Italy
 T. +39 0546 622063
 F. +39 0546 621640
 www.temasinergie.com

br


EXAMPLE 1

Relation between a 185cm wide corridor and the entrance through a 120cm wide door


EXAMPLE 2

Relation between a 140cm wide corridor and the encumbrance during the curve


EXAMPLE 3

Minimum useful size doors

External electric board **ONLY** for UL certification version - Dim. 25x60x100cm
For CE certification electric board is integrated

FOR ANY FURTHER QUESTIONS OR CONCERNS, PLEASE CONTACT TEMA SINERGIE

MAXIMUM WEIGHT TO BE MOVED: EUROPALLET (120x80cm) 1200 kg

Dimensions in [cm]

DESCRIPTION Manipulation Hot Cell Mod. FLEX 2T		TEMA SINERGIE Tema Sinergie S.p.A. Via Malpighi, 120 48018 Faenza (RA) Italy T. +39 0546 822963 F. +39 0546 821640 www.temasinergie.com	
DRAWING CODE SLST116 - 06	REV 00	STANDARD CONFIGURATION All pictures and specifications included in this document are indicative and may be subject to change without notice	
SCALE 1:50 - 1:25	DATE 04/10/2018	CUSTOMER	
DRAFTSMAN M.M.	APPROVED BY AS BUILT	All pictures and specifications included in this document are indicative and may be subject to change without notice	
DATE	SIGN	change without notice	
DRAFTSMAN	APPROVED	change without notice	
DATE	CHANGES	change without notice	
			FORMAT A3

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Faenza, September 11th 2018
Faenza, 11. september 2018

Your Re./ Vaša referencia		
Our Re./ Naša referencia	QUOTATION # 181205-2/A10 / CENOVÁ PONUKA č. 181205-2/A10	FG
SUBJECT/ PREDMET	New Hot Cell for BIONT – Bratislava / Nová horúca komora pre BIONT - Bratislava	

We hereby declare that the Hot cell FLEX 2T SAS fulfills the following requirements:

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ISO 14644-1, ISO 14644-7 a ISO 10648-2


We hereby declare also that Hot cell FLEX 2T SAS is GMP compliant.

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Find attached also a copy of the equipment declaration of conformity.

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**EC DECLARATION OF CONFORMITY OF THE
MACHINERY**

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